

Fast SeeBand

Single Step Protein Staining Solution

Warning: wearing gloves is highly recommended when handling SDS-PAGE and **Fast SeeBand** staining solution.

Storage: at room temperature.

This product is guaranteed for one year from date of purchase when properly handled and stored.

Sensitivity: 10 ng protein per band (BSA).

Important: This protocol is adjusted for 8 x 10 cm mini gel 1.4 mm thick.

Standard Protocol

Caution: Use caution while using the stain in a microwave oven. Do not overheat the staining solutions.

1. After electrophoresis, cut the stacking gel.
2. Put the gel in a container and add 25 ml of Fast SeeBand.
3. Shake the gel on an orbital shaker for 15 minutes.

Optionally: microwave on High Power for 30 seconds and shake 5 min.

Note: Stain intensity reaches a maximum within approximately 1 hour. Gels may be stained overnight without increased background.

4. Washing the gel with water can enhance staining sensitivity and increase staining over background.

Note: After using **Fast SeeBand** protein staining solution, close tightly the cap. Do not freeze.

Ordering information:

Cat. No.	Description
SB050	Fast SeeBand protein staining solution, 500 ml. Sufficient for 20 mini gels
SB051	Fast SeeBand protein staining solution, 60 ml. Sufficient for 2 mini gels