




## ClearBand ECL, Western Blotting Substrate

### Applications

- Western Blotting
- Dot Blot Assay

### Specifications

- Signal intensity: **Standard Grade**
- Compatible with Nitrocellulose & PVDF membrane
- For detection of horseradish peroxidase (HRP) enzyme activity



Detection with ClearBand ECL Western Blot substrate

#### ORDERING INFORMATION

| Catalog Number | Description   |
|----------------|---|
| CB-60          | Detection Reagents A and B, 30 ml each, good for <b>600 cm<sup>2</sup></b> membrane   |
| CB-120         | Detection Reagents A and B, 60 ml each, good for <b>1200 cm<sup>2</sup></b> membrane  |
| CB-250         | Detection Reagents A and B, 125 ml each, good for <b>2500 cm<sup>2</sup></b> membrane |

לפרטים ולהזמנות

[www.tivanbiotech.com](http://www.tivanbiotech.com)

עמית בקר: 054-2226019, פקס: 077-4701591

Email: [customerservice@tivanbiotech.com](mailto:customerservice@tivanbiotech.com)