


Blot&Go Hands Free, Autonomous Western Blot System

- From membrane to western in **less than 3 hours!**
- Can be left for running over night.
- Single step process: **just load the kit & go!**
- Support the use of Nitrocellulose and PVDF membranes.
- Single use, affordable and simple to use.
- Extended shelf life.
- Run 2 membrans at the same time.

no shaking • no electricity • no hassle

For orders:

CAT	Product
BNG-02	Blot&Go sample pack, 2 units, (including Go Buffer)
BNG-08	Blot&Go, 8 units (including Go Buffer)
BNG-16NB	Blot&Go, 16 units (without Go Buffer)
BNG-RL	Blot&Go, Roller
BNG-GB500	Go Buffer, 500 ml



Scan me with your smart phone to watch the video.
(or contact us at customerservice@tivanbiotech.com for direct link)



לפרטים ולהזמנות

www.tivanbiotech.com

עמית בקר: 054-2226019, פקס: 077-4701591

Email: customerservice@tivanbiotech.com